



# Ganj HAMG23PD Q61B

## 23 Litre B & S Class Autoclave



The 23L Ganj Autoclave is an entry-level, cost-effective compact steriliser that uses traditional dry vacuum. The HAMG23PD features imported high-quality parts including a **German** vacuum pump and an **Italian** water pump. Its **efficient condensing system** ensures optimal drying, while the **patented independent steam generator** guarantees thorough sterilisation.

The autoclave also features intelligent **load detection** for energy efficiency, a **large storage capacity** for up to **30,000** cycles, and comprehensive safety measures including **automatic water quality monitoring** and emergency stop functions. The HAMG23PD delivers reliable, **high-quality sterilisation** in a **user-friendly design**.



Fast Cycle Time From **28 Minutes**



**9 Programs:** Cycle Parameters Can Be Customised To Meet Various Needs



**User Friendly Design:** Intuitive design combined with a 3.5" LCD display.



**Large Storage Capacity:** Can Store Up To 30,000 Sterilisation Records



**Intelligent Load Detection:** Calculates & Uses The Least Amount of Resources

### Additional Features

- 3x Pre-Vacuum
- Pre-Heat Feature
- Vacuum Drying
- Intelligent Load Detection
- Embedded Control System
- Water Quality Monitor
- Intelligent Error Detection
- 3x Trays 395x190x23mm
- Built-In Large Memory
- Software Upgrade with USB
- Manual Lock
- Emergency Functions
- Power Off Protection

### High Quality Imported Parts

The series boasts imported high-quality parts, including a German Thomas vacuum pump (that can achieve a vacuum level of -0.90bar) and Italian solenoid valves and water pump, provide superior sterilization and maximum preservation for your instruments.

### Cycle Durations (1kg Load & Including Drying Time)

GEN134	Quick B	Quick S	user1	user2
49 min	34 min	24 min	as needed	as needed

Visit Our Website



### Optional Trays

The HAMG23PD comes default with a set of 3 trays. If you require more space to lay items on, an optional set of 5 trays can be purchased separately.

### Optional External Printer

An optional Hatmed cycle record external printer can be easily connected at any time. No service technician needed.

